

Name: Jiangxi Wellanimal Gene Technology Co., Ltd.
Address: 5F, Building 2, Ganzhou National High-level Talents
Science and Innovation Park(Phase II), 1 Wudang Mountain Road,
Zhanggong District High-tech zone, Ganzhou City, Jiangxi Province.

Telephone: 400-6723-898 E-mail: biolight\_info@163.com Website: http://www.wellanimal.com

## Anti-PCNA Antibody, Rabbit Polyclonal

Synonyms: Cyclin, DNA polymerase delta auxiliary protein, HGCN8729, MGC8367

Basic Info

Catalog

CP00140HuA10

Host

Rabbit

Con jugate

None

Size

100 μ L

Concentration

1.0 mg/mL

Physical State

Liquid

Species Reactivity

Human

Immunogen

Recombinant human PCNA protein,

fragment Met1~Ser261;

UniprotKB: P12004

Purification

Antigen Affinity Chromatography

Applications

WB/IHC/ICC

**Property** 

Form & Buffer: Supplied in PBS, 50% glycerol, PH7.4.

Specificity / Sensitivity: Anti-PCNA Antibody, Rabbit Polyclonal

recognizes endogenous levels of total PCNA protein.

Usage and

Storage

Shipped at  $4^{\circ}$ C.

Store at  $4^{\circ}C$  for frequent use.

Aliquot and store at -20℃ for 12 months.

Avoid repeated freezing/thawing and violent shaking.

Please centrifuge it, before using.

**Applications** 

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


WB:  $0.5 \sim 5 \mu \text{ g/mL}$ 

IHC: 5~20 μ g/mL

IF: 5~20 μ g/mL

ICC:  $5\sim20 \mu \text{ g/mL}$ 

QC Data




## Figure 1. Application in WB

Western blot analysis of extracts of various cell lines, using PCNA antibody (CP00140HuA10) at  $0.5\,\mu\,g/mL$ . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Blocking buffer: 5% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 7 s.



QC Data



## Figure 2. Application in IHC

Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) analysis of Canine breast cancer tissue sections labelling PCNA with purified CP00140HuA10 at 5  $\mu \text{g/mL}.$  Heat mediated antigen retrieval was performed using citrate buffer (pH 6.0). Tissue was counterstained with Hematoxylin. Rabbit specific IHC polymer detection kit HRP/DAB secondary antibody was used at 1/4000 dilution. PBS instead of the primary antibody was used as the negative control.

For Research Use Only!